

FIGGERT
ENDOTECH



URO



 **MADE IN GERMANY**



ACC



LAP



GYN




URO



ART



ORL



REP

The following brochures are available:

- *Devices + Accessories*
- *Laparoscopy*
- *Gynecology*
- *Urology*
- *Arthroscopy*
- *ENT (ORL)*
- *Rep.-Service Endoscopes*

Certificates see last page



URO - Urology

<u>Table of Contents</u>		Page
Index		1
<hr/>		
<u>Endoscopes:</u>		
	Ø 4 mm - <i>HD-Line</i> & <i>A-Line</i> (Standard) & <i>X-Line</i> (Low-Cost)	2
	Ø 2,0 mm - 2,9 mm & other connections (Ø 4 mm)	3
<hr/>		
<u>Resectoscopy:</u>		
	Working-Elements for Teens & Adults	4
	Resectoscope Sheaths for Teens & Adults	5
	Monopolar Electrodes - Single & Double Stem	6
	Resectoscope for Children incl. Cystoscope & Urethrotome	7
<hr/>		
<u>Uretero-Renoscopy:</u>		
	Uretero-Renoscopes - Semi-Rigid	8
	Instruments & Electrodes for Uretero-Renoscopes	9
<hr/>		
<u>Cystoscopy:</u>		
	Cysto-Urethroscopy Sheaths CH / FR 7,5 - 25	10
	Endoscope Bridges - Albarrans with 1 & 2 channels	11
	Flexible Instruments - WL 40 cm	12
	Electrodes for Coagulation - Monopolar	13
	Optical Instruments for Cystoscopes with Ø 4 mm	14
	Urethrotomy - Sheath CH / FR 21	15
	Accessories - Bladder Syringe & Ellik-Evacuator	16
	Accessories - Stone Extractor & Polypectomy Snare	17
<hr/>		
Nephroscopy	Percutaneous Nephroscopy	18
Lithotripsy	Mechanical Lithotripter	19

Endoscopes for Urology - Ø 4,0 mm



Cystoscopes - Hysteroscopes - **HD-Line** - HIGH DEFINITION:

Ø	Working length	Direction of view	REF.-NO. :
4,0 mm	30 cm	0°	300.0400.HD
		12°	300.0412.HD
		30°	300.0430.HD
		70°	300.0470.HD

Cystoscopes - Hysteroscopes - **A-Line** (Standard):

Ø	Working length	Direction of view	REF.-NO. :
4,0 mm	30 cm	0°	300.0400.A
		15°	300.0415.A
		30°	300.0430.A
		70°	300.0470.A

Cystoscopes - Hysteroscopes - **X-Line** (Low-Cost):

Ø	Working length	Direction of view	REF.-NO. :
4,0 mm	30 cm	0°	300.0400.X
		30°	300.0430.X
		70°	300.0470.X

Endoscopes for Urology - Ø 2,0 mm - 2,7 mm - 2,9 mm



Cystoscopes - Hysteroscopes - **A-Line** (Standard):

Ø	Working length	Direction of view	REF.-NO. :
2,0 mm	26 cm	0°	200.0200.26
		30°	200.0230.26
	30 cm	0°	200.0200.A
		30°	200.0230.A
2,7 mm	30 cm	0°	301.2700.A
		30°	301.2730.A
		70°	301.2770.A
2,9 mm	30 cm	0°	301.2900.A
		30°	301.2930.A
		70°	301.2970.A

Cystoscopes - Hysteroscopes - **A-Line** - Richard Wolf ® connection :

Ø	Working length	Direction of view	REF.-NO. :
4,0 mm	29 cm	0°	300.0405.A
		12°	300.0412.A
		30°	300.0425.A
		70°	300.0475.A

Cystoscopes - Hysteroscopes - **A-Line** - Olympus ® connection :

Ø	Working length	Direction of view	REF.-NO. :
4,0 mm	28,0 cm	0°	307.0400
		12°	307.0412
	28,2 cm	30°	307.0430
		28,3 cm	70°

Working Elements - for Teens & Adults

TURP - Trans-Urethral Resection of the Prostate:

Working element " ACTIVE " :


for	FR	f. Scope	REF.-NO. :
Teens	19/22	2,7 / 302	350.0021
Adults	24/26	4,0 / 302	350.0011
Children	see page 7		



Urethrotomy - Hystero-Resection:

Working element " PASSIVE " :

for	FR	f. Scope	REF.-NO. :
Teens	19/22	2,7 / 302	350.0020.1
Adults	24/26	4,0 / 302	350.0010.1.OP
Children	see page 7		



Electrosurgical Cables - Monopolar - 3 m

Plug of cable for Electrode



Plug of cable for Electrosurgical Units of...

REF.-NO. :



...EMED ®, Aesculap ®, Martin ®,
Siemens ®, Erbe ® (T-Models)

4mm plug-Ø

350.0001.N



...Karl Storz ®, Richard Wolf ®,
Erbe ® (ICC/ACC-Models)

5mm plug-Ø

350.0001.S




...Valleylab ®, Bowie ®,
Eschmann ®, Keymed ®

8 mm plug-Ø

350.0001.V

Resectoscope Sheaths - for Teens & Adults



Continuous-Flow Sheath - complete with Obturator

for	FR	Endoscope	not rotatable	rotatable :
Children	13	2,7 / 187,5	see page 7	REF.-NO.:
Teens	19/22	2,7 / 302	360.1922	360.1922.R
Adults	24/26	4,0 / 302	360.0426	360.2426
High-Flow	27/28	4,0 / 302	360.0428	360.2428



Deflecting Obturator

CH/FR 24/26	REF.-NO.:	360.0824
CH/FR 27/28	REF.-NO.:	360.0927











Visual Obturator

	<u>without channel</u>	<u>with channel</u>
CH/FR 24/26	360.0823	360.0823.K
CH/FR 27/28	360.0830	360.0830.K

Monopolar Electrodes - Single & Double Stem



		CH / FR 19 / 22 Single Stem	CH / FR 24 / 26 Single Stem	CH / FR 24 / 26 Double Stem
	Loop - retrograde curved	REF.-NO.: 355.1924	REF.-NO.: 352.0924	REF.-NO.: 353.0924
	Loop - straight	355.1910	352.0910	353.0910
	Ball	3 mm 5 mm 355.1903 -	352.0902 352.0900	353.0902 353.0900
	Roller barrel	3 mm 5 mm 355.1906 -	352.0906 352.0916	353.0906 353.0916
	Roller barrel serrated	3 mm 5 mm -	352.0930 352.0932	353.0930 353.0932
	Cone	-	352.0904	353.0904
	Knife	355.1912	352.0912	353.0912
	Needle	-	352.0918	-


Resectoscope for Children - incl. Cystoscope & Urethrotome

Resectoscope, one fixed stop cock,
not rotatable, consisting of:

Working element
Resectoscope sheath
Obturator

REF.-NO. :

360.2713.AE
360.2713.S
360.2713.O



Electrodes:

Ball	Ø 2 mm		360.2730
Knife	90°		360.2731
Loop	90°	small	360.2732
		large	360.2733
Hook			360.2734

360.2730
360.2731
360.2732
360.2733
360.2734



Cystoscope Sheath with working channel 3 FR
WL 136 mm, inclusive Obuturator
(Instruments FR 3 on request)

360.2720.S



Telescope bridge to cystoscope sheath

360.2721



Albarran catheter deflecting mechanism

360.2722



Urethrotome sheath FR 13
Obturator
Telescope bridge

360.2740
360.2740.O
360.2741



Cold knives	straight	360.2742
	round	360.2743
	hook	360.2742
	Sickle	360.2743

360.2742
360.2743
360.2742
360.2743



Flexible electrodes FR 3	button	330.0301
	needle	330.0302

330.0301
330.0302



for endoscopes Ø 2,7 mm - working length 187,5 mm:

Direction of view	0°	502.2700.A
	30°	502.2730.A

502.2700.A
502.2730.A



Uretero-Renoscopes - Semi-Rigid

Features:

- CAD-designed tension-free composed-objective for optimal sharpness, colour sync and resolution !
- Image transmission through optical fiber-bundle for high resolution and optimal colour sync !
- High quality light fibers !
- Optical glass: Distal and proximal Sapphire !
- Direction of view 12° (to working channel)
- Autoclavable
- Two stop cocks for suction and irrigation
- Two insert possibilities for instruments



Standard-Line - with up to 30.000 Pixel

Two different eyepiece versions available: lateral (.B) and laterally offset (.M)

Size	Ø - Steps from distal tip to proximal end		Working Channel		Working Length	
	Ø in mm	(FR)	Ø	(FR)	315 mm	430 mm :
S	2,0 - 2,5 mm	(6,0 - 7,5)	1,4 mm	(4,2)	365.6031.M 365.6031.B	365.6043.M 365.6043.B
M	2,7 - 3,3 mm	(8,0 - 9,8)	1,8 mm	(5,4)	365.8031.M 365.8031.B	365.8043.M 365.8043.B
L	2,8 - 3,8 mm	(8,5 - 11,5)	2,1 mm	(6,3)	365.8531.M 365.8531.B	365.8543.M 365.8543.B



HD-Line - with up to 40.000 Pixel

Size	Ø - Steps from distal tip to proximal end		Working Channel		Working Length		
	Ø in mm	(FR)	Ø	(FR)	150 mm	310 mm	425 mm :
S	1,6 - 2,7 mm	(4,9 - 8,1)	1,33 mm	(4,0)	365.0127	365.0327	365.0427
M	2,0 - 3,0 mm	(6,0 - 9,0)	1,70 mm	(5,1)	365.0130	365.0330	365.0430
L	2,5 - 4,0 mm	(7,5 - 12,0)	2,26 mm	(6,8)	365.0140	365.0340	365.0440
XL	2,5 - 4,4 mm	(7,5 - 13,2)	2,26 mm	(6,8)	-	-	365.0444

Sterilization basket for URS - Dimensions 640 x 150 x 77 mm










961.6415

Instruments & Electrodes - for Uretero-Renoscopes






Semi-Rigid Instruments - WL 60 cm

	CH / FR	3	4	5
		REF.-NO.:	REF.-NO.:	REF.-NO.:
	Biopsy forceps	366.0301	366.0401	366.0501
	Foreign body forceps	-	366.0403	366.0503
	Stone street	-	366.0404	366.0504
	Ureter strictures scissors	-	366.0405	366.0505
	Alligator forceps	366.0307	366.0407	366.0507
	Large fragments	-	366.0408	366.0508
	3-prong retriever, flexible	365.0317	-	365.0517

Electrodes for Coagulation - Monopolar - WL 65 cm



		CH/FR 5	CH/FR 7
		REF.-NO.:	REF.-NO.:
	Ball	330.0501.L	330.0701.L
	Needle	330.0502.L	330.0702.L
	Loop	330.0504.L	330.0704.L

Cysto-Urethroscopy Sheaths

Sheath with
Standard Obturator



Visual Obturator



Ø CH / FR	Colour Code	Sheath with Standard Obturator		Compatible Endoscope
		REF.-NO. :		REF.-NO. :
Newborns, Infants, Children:			Working Channel	
7,5	Yellow	360.1975	-	501.1930
9,5	Blue	360.2095	CH / FR 3	200.0230.187
11	Green	360.2110	CH / FR 3	200.0230.187
13	Red	360.2720.S	CH / FR 3	502.2730.A (Bridge 360.2721)
15	Black	310.0015	CH / FR 5	301.2930.A
Adults:			Visual Obturator	
17	Yellow	310.0017	311.0017	300.0430.A
19	Green	310.0019	311.0019	300.0430.A
21	Red	310.0021	311.0021	300.0430.A
23	Blue	310.0023	311.0023	300.0430.A
25	White	310.0025	311.0025	300.0430.A

Endoscope Bridges - Albarrans

- 310.0100** Endoscope bridge for diagnostic
- 310.0101** Endoscope bridge with 1 lockable instrument channel
- 310.0102** Endoscope bridge with 2 lockable instrument channels



Endoscope bridges	Cysto-Urethroscope-Sheath sizes in FR/CH (see page 13)				
	Sheath FR 17	Sheath FR 19	Sheath FR 21	Sheath FR 23	Sheath FR 25
REF.-NO.	for use with flexible instruments of following diameters:				
310.0101	1 x 5 FR	1 x 6 FR	1 x 7 FR	1 x 9 FR	1 x 10 FR
310.0102	- 1 x 5 FR	2 x 4 FR 1 x 6 FR	2 x 5 FR 1 x 7 FR	2 x 6 FR 1 x 9 FR	2 x 8 FR 1 x 10 FR

- 310.0111**
Albarran catheter deflecting mechanism, 1 instrument channel, lockable






- 310.0112**
Albarran catheter deflecting mechanism, 2 instrument channels, lockable



Albarran deflecting mechanism	Cysto-Urethroscope-Sheath sizes in FR/CH (see page 13)				
	-	Sheath FR 19	Sheath FR 21	Sheath FR 23	Sheath FR 25
REF.-NO.	for use with flexible instruments of following diameters:				
310.0111	-	1 x 5 FR	1 x 6 FR	1 x 9 FR	1 x 10 FR
310.0112	-	2 x 4 FR 1 x 5 FR	2 x 5 FR 1 x 6 FR	2 x 6 FR 1 x 9 FR	2 x 8 FR 1 x 10 FR




Flexible Instruments - WL 40 cm



		5 CH/FR	7 CH/FR
		REF.-NO.:	REF.-NO.:
	Biopsy forceps	320.0501	320.0701
	Scissors, pointed, single action	320.0502	320.0702
	Foreign body forceps	320.0503	320.0703
	Biopsy forceps with 3 spikes	320.0504	-
	Scissors, blunt, single action	320.0505	320.0705
	Hook scissors, single action	320.0506	320.0706
	Alligator grasping forceps	320.0507	320.0707
	Punch, through-cutting	320.0508	-

Electrodes for Coagulation - Monopolar



		<u>WL</u>	<u>CH/FR 5</u>	<u>CH/FR 7</u>
			<u>REF.-NO.:</u>	<u>REF.-NO.:</u>
	Ball	45 cm	330.0501	330.0701
		65 cm	330.0501.L	330.0701.L
	Needle	45 cm	330.0502	330.0702
		65 cm	330.0502.L	330.0702.L
	Loop	45 cm	330.0504	330.0704
		65 cm	330.0504.L	330.0704.L

Electrosurgical Cables - Monopolar - 3 m



Plug of cable for Electrode



Plug of cable for Electrosurgical Units of...

REF.-NO. :



...EMED ®, Aesculap ®, Martin ®,
Siemens ®, Erbe ® (T-Models)

4mm plug-Ø

350.0001.N



...Karl Storz ®, Richard Wolf ®,
Erbe ® (ICC/ACC-Models)

5mm plug-Ø

350.0001.S



...Valleylab ®, Bowie ®,
Eschmann ®, Keymed ®

8 mm plug-Ø

350.0001.V

Optical Instruments - for Cystoscopes with Ø 4 mm



for Cystoscopes Ø 4 mm - 0° + 30° and
Cysto-Urethroscope-Sheaths FR 21, 23, 25 :

REF.-NO.:



Biopsy forceps

312.0011



Biopsy punch

312.0012



Scissors

312.0013



Grasping forceps

312.0014



Stone grasping forceps
single action

312.0054

for Cystoscopes Ø 4 mm - 30° + 70° :







REF.-NO.:



Biopsy forceps

312.0015

Urethrotomy

	REF.-NO.:
Urethrotome-Sheath CH / FR 21 with instrument channel CH / FR 5	355.0021
	
Obturator to Urethrotome-Sheath	355.0026
	
Guide tube to urethr.-sheath, lateral open for inserting catheters	355.0022
Continuous-Flow sheath for urethrotom CH / FR 21	355.0023
Telescope bridge to urethr.-sheath, Channel CH / FR 5	355.0031
	
Cold knife, straight	352.1011
	
Cold knife, round	352.1012
	


Urethrotome acc. to Otis-Mauermayer, complete in wooden case, consisting of:
3 tips: conical, flat, ball; spare knife

357.0011 Spare knife:



357.0001

Bladder Syringe & Ellik-Evacuator



Bladder syringe acc. to Reiner-Alexander	REF.-NO.:	Spare barrel
50 ml	315.6005	315.6006
75 ml	315.6007	315.6008
100 ml	315.6010	315.6011
150 ml	315.6015	315.6016
Conical adapter for catheters	315.6002	

Ellik-Evacuator - complete
with rubber ball and glass box

SET-NO. :

Set with plug cone

315.5001

Set with locking cone

315.5011



Spare parts

REF.-NO. :

Rubber ball

315.5003

Glass box

315.5002

Plug cone

310.0104

Locking cone

310.0105

Stone Extractor & Polypectomy Snare



Handle without HF-Connection

CH / FR 3 + 5
CH / FR 7

REF.-NO.:

330.1600
330.1600.7



Stone Extractor - 4 wires - WL 60 cm - without handle

CH / FR 3
CH / FR 5
CH / FR 7

330.0312
330.0512
330.0712

Stone Extractor - 6 wires - WL 60 cm - without handle

CH / FR 5
CH / FR 7

330.0512.6
330.0712.6



Handle for Polypectomy Snare, with LL-adapter and HF-connection

REF.-NO.:

330.0600



Polypectomy Snare - WL 60 cm - without handle

CH / FR 5
CH / FR 7

330.0510
330.0710

Percutaneous Nephroscopy (Extracts only - Complete Set on request)

REF.-NO. :



305.0406.A

Nephroscope with parallel eyepiece
Direction of view 0°, autoclavable
working channel Ø 5 mm, with insertion tube,
add. two flanking working channels
with closable stopcocks for 1mm-Instruments



OP-Sheath, 1 fixed stop cock

CH/FR 24

380.0024.A

Hollow-Obturator to OP-Sheath

380.0024.O



OP-Sheath, 1 rotatable stopcock

CH/FR 26

381.0026

Hollow-Obturator to OP-Sheath

381.0026.O



Instruments, semi-rigid WL 40 cm

9 CH/FR:



Biopsy forceps

380.0901



Foreign body forceps

380.0903



Scissors

380.0905



Hook scissors

380.0906.



3-Arm-Grasper
(U-Handle) - WL 33 cm

380.0917.9CH

Compatible to sheath 381.0026 only:

14 CH/FR:



Rigid 3-Arm-Grasper (U-Handle)
for large stones - WL 33 cm

380.0917

Lithotripsy - Mechanical



REF.-NO. 313.2301

Lithotripter Punch Ø CH / FR 23

2 stop cocks with LL-adapter

for Cystoscope with direction of view 30°: REF.-NO. 300.0430.A



REF.-NO. 313.2401

Lithotripter Crusher Ø CH / FR 24

mechanical, with handle and turning srew, 2 stop cocks with LL asapter,

for Cystoscope with direction of view 70°: REF.-NO. 300.0470.A

URO



ACC



LAP



GYN



ART



ORL



REP



DEU:

*Fiegert-Endotech Deutschland
78532 Tuttlingen, Germany
fon: (00 49) 74 62 / 94 85 - 0
fax: (00 49) 74 62 / 94 85 - 22
info@fiegert-endotech.com
www.fiegert-endotech.com*

USA:

*Fiegert-Endotech-Inc. U.S.A.
Ft. Lauderdale, Florida 33312
fon: (00 1) 954 748 2700
fax: (00 1) 954 748 2743*